

## OHAUS SOLUTION BALANCE

The heavy duty solution balance is the industry standard for precision in high-capacity weighing. The front beam has a weighing capacity of 100 g while the back beam has a capacity of 1000 g. The locking sliding tare weight compensates for containers up to 2270 grams. Included are auxiliary weights for weighing over 1100 g. The extra weights are store on rack when not in use. Also features ; magnetic dampening, front and rear parallax free reading device, notched beams, precision ground knives and an 11" stainless steel weighing plate.

MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	TARE
G4421	HEAVY DUTY SOLUTION BALANCE	20 KG	1 G	11" DIA.	2270 G



ACCESSORIES	DESCRIPTION
G4422	SCOOP & WEIGHT
G4424	EXTRA WEIGHT, 1 KG
G4425	EXTRA WEIGHT, 2 KG
G4426	EXTRA WEIGHT, 5 KG
G4427	EXTRA WEIGHT, 10 KG
G4428	EXTRA WEIGHT SET

## OHAUS FIELD SCALE

The Ohaus Field Scale is designed is an affordable compact scale for laboratory, field, and general usage. The large, stabilized pan can hold oversized samples without tipping, and the weighted base keeps the entire scale stable. Reading are taken from a single, sliding-weight beam, while a set of slotted weights, stored in the base can extend the capacity. Available in English and Metric standards. The scale has a built in carrying handle and 11" platform.

MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM
G5111	FIELD SCALE	36 LBS	0.01 LBS	11" DIA.
G5121		16 KG	5 G	



ACCESSORIES	DESCRIPTION
G5112	HANGER
G5115	EXTRA WEIGHT SET, (ENGLISH)
G5125	EXTRA WEIGHT SET, (METRIC)
G5117	EXTRA WEIGHT, 1 LB
G5118	EXTRA WEIGHT 2 LBS
G5119	EXTRA WEIGHT 5 LBS
G5120	EXTRA WEIGHT 10 LBS