

OHAUS HANDHELD SCALES



These precision Ohaus Hand-Held Scales are convenient and compact allowing weighing in the palm of your hand. Constructed from rugged plastic these scales weigh in g, oz, dwt, and oz t. Some features include an easy to read LCD display and an auto calibration achieved from the keypad. Powered by 2 AA alkaline batteries (included) or by an optional AC adapter. A built-in energy saving auto shut-off senses inactivity after two minutes. Perfect in the laboratory and in the field these scales are ideal to carry anywhere for quick and accurate readings.

MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	UNITS
G4500	HANDHELD SCALE	120 G	0.1 G	3.25" X 3.0"	G, OZ, DWT, OZ T
G4502		120 G / 60 G	0.2 G / 0.1 G		
G4504		320 G	0.1 G		

OHAUS COMPACT SCALES

The Ohaus Compact Series Scales offer accurate measurements at an economical price. The sleek, low profile design constructed from rugged plastic makes these scales the perfect portable choice in the lab or the field. Features include a large high contrast display with an energy-saving auto shut-off after 5 minutes of inactivity (battery powered mode only). Can operate from 3 AA alkaline batteries (included) or optional AC adapter. Scale stabilizes in 3 seconds while digital calibration is achieved from the keypad.



MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	UNITS
G4520	COMPACT SCALE	200 G	0.1 G	5.75" X 5.25"	G, LBS: OZ
G4522		2000 G	1 G		
G4524		5000 G	2 G		

OHAUS SCOUT II SERIES BALANCE

The Scout is an affordable alternative with superior performance available in 200 to 6000 gram capacities. The series has a simple 2-button operation with a large, high contrast LCD display with a 3 second stabilization time. Also featured, is a large stainless steel weighing surface and weigh below hook. Multiple weight units are user selectable- g, oz, oz t, dwt, and parts counting. Calibration is made easy digitally from the keypad. The balance is powered by either 100, 120, 220, 240V AC, 50/60 Hz adapter which is furnished or one 9V battery (not included). Auto shut-off and low battery indicator are standard.



MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	UNITS
G4525	SCOUT II BALANCE	200 G	0.01 G	6.75" X 2.4" X 8.25"	G, LBS : OZ
G4528		400 G	0.1 G		
G4530		600 G			
G4532		1200 G			
G4534		6000 G	1 G		
G4536		2000 G			