

WEIGHTS AND WEIGHT SETS

The metric weight sets are used for calibration. Weights are constructed of stainless/aluminum and are adjusted to NIST (NBS) Class C Tolerance. Hogentogler will issue a Certificate of Traceability upon request. The cost of the certificate is based on the number and tolerance of weights certified. Weight sets are furnished with cases.



MODELS	DESCRIPTION	STANDARDS
G5511	WEIGHT SET, 100 G X 1 MG	NIST (NBS) CLASS C TOLERANCE
G5661	WEIGHT SET, 500 G X 1 G	
G5662	WEIGHT SET, 1,000 G X 1 G	
G5663	WEIGHT SET, 2,000 G X 1 G	
G5663T	WEIGHT SET, 2,000 G X 1 G, W/ CERTIFICATE	

ACCESSORIES	DESCRIPTION	SIZE
G4681	CLASS 1 CALIBRATION WEIGHT	100 G
G4682		200 G
G4683		500 G
G4684		1,000 G
G4685		2,000 G
G4688		4 KG
G4689	5 KG	
G4686	CLASS 4 CALIBRATION WEIGHT SET	500 G X 100 MG
G4687		2 KG X 100 G

TEST WEIGHTS, CLASS F

These Class F Test Weights are for checking the accuracy of scales. The weights are slotted with a handle and are constructed from cast iron.



MODELS	DESCRIPTION	STANDARDS
G4691	TEST WEIGHT, 5 LB	CLASS F
G4692	TEST WEIGHT, 10 LB	
G4693	TEST WEIGHT, 25 LB	
G4694	TEST WEIGHT, 50 LB	
G4695	TEST WEIGHT, 100 LB	
G4701	TEST WEIGHT, 5 KG	
G4702	TEST WEIGHT, 10 KG	
G4703	TEST WEIGHT, 25 KG	
G4705	TEST WEIGHT, 50 KG	
G4706	TEST WEIGHT, 100 KG	

UNIT WEIGHT SCALES

These full capacity portable bench beam scales are for field or lab tests of density and unit weight. These self-contained scales are compact, portable, easy to operate, and rugged for use over extended periods of time. The scales housing is constructed of all steel with floating platform for lasting accuracy. Special cast aluminum alloy levers with adjustable hardened steel pivots. Self-aligning hardened steel bearings with friction points. These scales are furnished with mirrored balance indicators for fast, accurate and convenient viewing. Also, the reading beam is notched with clear and sharp graduations. The scales have a built-in platform and beam locking mechanisms to prevent damage during transport or during idle periods. They also feature a bubble level indicator and adjustable leveling legs. Each scale is factory calibrated and subject to rigid quality control standards. A carrying handle is also available (**G5376** or **G5377**).

**G5371-G5374****G5373 & G5375**

MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM
G5371	UNIT WEIGHT SCALES	52 LBS	.01 LB	12" X 17"
G5372		131 LBS	.01 LB	
G5373		301 LBS	1 LB	
G5374		26 KG	.01 KG	
G5375		52 KG	.02 KG	

ACCESSORIES	DESCRIPTION
G5376	CARRYING HANDLE FOR G5371, 72, 74 & 75
G5377	CARRYING HANDLE FOR G5373