

OHAUS NAVIGATOR SERIES BALANCE

These portable high precision balances are ideal for almost any application including parts counting. They have fifteen standard weighing units including grams, ounces, ounces troy, carats, and pounds, as well as one custom weighing unit. A high contrast LCD displays weight and text messages and easy to use interface with cursor keys to navigate menus and functions are just some of the outstanding features. Each balance is furnished with an AC adapter or can use 8 AA Alkaline batteries (not included). A bi-directional RS232 interface is also available as an accessory to utilize the advanced options. Other features include, auto zero tracking, stability levels, auto calibration, high speed display update, and auto shut-off.



MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	UNITS
G4596	NAVIGATOR SERIES BALANCE	32 G	0.002 G	10.5" X 7.25" X 2.25"	G, MG, KG, OZ T N, OZ, TALES (3) DWT, CT, TI, M TOLA, CUSTOM
G4598		210 G			
G4600		410 G	0.01 G		
G4602		810 G			
G4603		2100 G			
G4605		4100 G	0.1 G		
G4607		8100 G			

ACCESSORIES	DESCRIPTION
G4610	RS232 BI-DIRECTIONAL INTERFACE

OHAUS ADVENTURER SERIES BALANCE

The Ohaus Adventurer Series Balance incorporates advanced construction and rugged design to deliver added durability. These top loading balances are designed with an impressive range of features which include, a simple three-button operation and ergonomic leveling system. Multiple weighing units including parts counting. RS232 interface is standard on all models while the **G4640** & **G4641** includes a draft shield.



MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	UNITS
G4640	ADVENTURER SERIES BALANCE	150 G	0.001 G	8.5" X 8.3" X 13.5"	G, KG, MG, CT N, LB, OZ, OZ T GN, DWT, MOMMES MESGHALS, TALES (3)
G4641		310 G			
G4642		510 G	0.01 G	3.9" DIA.	
G4643		1500 G			
G4644		3100 G			
G4645		4100 G	0.1 G	7.1" DIA.	

OHAUS CHAMP MULTI-FUNCTIONAL BENCH SCALES

When custom performance is needed, the Champ's multi-functional models offer versatility like parts counting, dynamic weighing percent, and checkweighing. These scales features a large 2-line LCD display with 8 function, 4 cursor membrane keys. The stainless steel platter with painted steel carbon frame. Bi-directional RS232 with GMP data output is also standard. Weigh modes include weighing, parts counting, dynamic weighing, percent weighing, and check weighing. Powered by AC adapter (supplied)



MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	UNITS
G4650	CHAMP MULTI-FUNCTIONAL BENCH SCALE	10 KG	0.0005 G	12" X 12" X 4"	LB, KG, G, CUSTOM UNIT
G4652		25 KG	0.001 G		
G4654		50 KG	0.002 G	18" X 18" X 4.1"	
G4656		100 KG	0.005 G		
G4658		250 KG	0.01 G		

OHAUS CHAMP GENERAL PURPOSE BENCH SCALES

The Ohaus Champ General Purpose Bench Scales offer custom versatility built on a rugged steel frame with leveling feet and leveling bubble. Equipped with a precise, easy-to-use indicator, featuring weighing and parts counting. A RS232 interface is standard on all models. Powered by the supplied AC adapter or 6 Alkaline "C" batteries (not included).



MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	UNITS
G4661	CHAMP GENERAL PURPOSE BENCH SCALE	10 KG	0.001 G	12" X 12" X 4"	LB, KG, G
G4663		25 KG	0.002 G	18" X 18" X 4.1"	
G4665		50 KG	0.005 G		
G4667A		100 KG	0.005 G		
G4667		100 KG	0.01 G		
G4669		250 KG	0.02 G		

OHAUS ES SERIES BENCH SCALES

These multipurpose electronic bench scales are economical and durable. Easy, two-button operation and large display guarantee accurate results. The TRIAX indicator stand allows you to mount the indicator to the scale base, above the scale platform, or on any vertical surface. These scales have large stainless steel platforms with a black painted frame. Powered by 9V DC, 500 mA AC adapter (included) or 6 Alkaline "AA" batteries (not included). The **G4670** carrying case is also available.



MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	UNITS
G4671	ES SERIES ALL-PURPOSE BENCH SCALE	30 KG	0.01 G	12.2" X 10.8" X 1.5"	KG, LB, OZ, LB:OZ
G4673		50 KG	0.02 G		
G4674		100 KG	0.05 G	20.5" X 15.7" X 3"	
G4675		200 KG	0.1 G		

OHAUS EXPLORER SERIES BALANCES

The Ohaus Explorer Series Balances have incredible combination of performance and value. Four cursor keys provide easy menu navigation and lead you through on-screen prompts for a wide range of applications. Also featured is a custom LCD display with seven alphanumeric characters. Easy to upgrade when additional features are needed. The **G4783** & **G4785** are supplied with a draft shield.



MODELS	DESCRIPTION	CAPACITY	READABILITY	PLATFORM	UNITS
G4783	EXPLORER SERIES BALANCE	210 G	0.001 G	4.7" DIA.	G, MG, KG, LBS, OZ, OZ T, CT, DWT, TAELS (3) MOMMES, GN TL, N, CUSTOM
G4785		410 G			
G4787		100 G / 410 G			
G4789		610 G	0.01 G	6.8" X 6.8"	
G4791		2,100 G			
G4793		4,100 G			
G4795		4,100 G	0.1 G	8" X 8"	
G4797		1000 G / 4,100 G	0.01 G / 0.1 G	6.8" X 6.8"	
G4802		6,100 G	0.1 G	8" X 8"	
G4635A		8,100 G			
G4803		12,000 G		11" X 14"	
G4804		22,000 G			
G4805		32,000 G			



HOGENTOGLER & CO., INC.

TOLL FREE 800-638-8582

www.hogentogler.com