

Customer Rebates!

July 1, 2022 – Nov. 30, 2022

seca
Precision for health

#1 Physician's Choice Digital Column Scale

\$25.00
customer
rebate!




seca 777 | Item # 7771821004

Digital Column Scale with eye-level display

- + **8-year warranty (4x industry standard!)**
- + Over-sized, low-profile full metal platform (no plastic)
- + Tip-resistant base
- + Eye-level display swappable to left or right side
- + Robust rubber coated castors for easy transport
- + Available in lbs/kg or kg-only
- + USB interface for PC connectivity

Capacity	550 lbs / 250 kg
Graduation	0.2 lbs / 100 g
Dimensions (WxHxD)	16.1 x 53.4 x 22.7 inch / 41.0 x 135.6 x 57.6 cm
Platform (WxHxD)	13.8 x 2.3 x 16.5 inch / 35.0 x 5.9 x 41.8 cm
Measuring range	4 – 90 inch / 10 – 230 cm
Weight	28.0 lbs / 12.7 kg
Power supply	Batteries, power adapter (optional)
Data transmission	USB interface for PC connectivity
Functions	TARE, HOLD, lbs/kg switch-over lock, BMI
Electromagnetic radiation resistance	60601-1-2 compliant
Accessories	+ Handrail seca 477 + Panda seca 459 + Power adapter 401

The king of stability!



\$10.00
customer
rebate!

seca 477 | Item # 4770000009

Optional Integrated Handrail

- + Enhanced patient safety without affecting measurement
- + Extra support and stability during the weighing process
- + Enables a firm stance for patients in need of assistance

Dimensions (WxHxD)	16.9 x 44.2 x 12.2 inch / 42.9 x 112.2 x 30.9 cm
Weight	8.8 lbs / 4.0 kg



Optional handrail for increased patient stability

To accommodate patients in need of extra support during the weighing process, the seca 777 is equipped with the option to add an integrated handrail (seca 477). The easy-to-mount handrail raises the level of patient security to a new level. The measuring rod is securely attached to the column and provides a wide measuring range that is suitable for pediatrics to tall individuals. In addition, we've added extra space between the column and the platform to ensure the highest comfort.



Cost-effective without any compromises in quality

The seca 777 puts unstable and inaccurate weighing and measuring to an end. It boasts a tilt-proof full metal platform and a stable single piece column. The eye-level display and ultra sturdy base provide comfort and security. It is a quality every patient will notice, as soon they sets foot on the platform. The scale is easy to move by staff due to the large rubber coated castors and a wide wheelbase.

To redeem your rebate, contact seca at info.us@seca.com
Proof of purchase must be included. Rebate must be redeemed within 30 days of purchase.

Precision for health
seca.com/777