# LRB23W1HC

# Mediline







#### Laboratory Refrigerator with Bottom Mount Compressor and Electronic Controls

LED



## **Unit Dimensions**

| 80 15/32"<br>2044 mm<br>2044 mm | 313/hr 806 mm) |
|---------------------------------|----------------|
| 4                               | ,              |

### Laboratory refrigerator

with electronic controls

Total gross/net capacity I (cu.ft.)

Exterior dimensions in mm (w/d/h)
in inches (w/d/h)

 $\begin{array}{c} \textbf{Interior dimensions} & \text{in mm (w/d/h)} \\ & \text{in inches (w/d/h)} \end{array}$ 

Energy consumption in 365 days \* Ambient temperature in °C (°F)

Refrigerant
Sound rating \*\*
Connection rating
Frequency/voltage

Cooling system
Defrost
Temperature range in °C (°F)

External cabinet finish

Door/lid material Internal liner material

Type of control
Temperature display
Alarm

Alarm forwarding
Interior light
Adjustable shelves
Shelf material

Shelf load weight in kg (lbs)
Usable shelf area in mm (inches)

Casters

Handle
Lock
Self-closing door
Door hinges
Insulation in mm (inches)
Gross/net weight in kg (lbs)

Accessories

Shelf, plastic-coated

#### LRB23W1HC

Mediline

664/627 (23.3/22.0) 747/750/2044 29 7/16"/29 9/16"/80 1/2" 650/581/1710 25 5/8"/22 7/8"/67 3/8"

344 kWh

10°C to 35°C (50°F to 95°F)

R290 55 dB(A) 3.0 A 60 Hz/115V~ forced-air

automatic 3°C to 16°C (37°F to 61°F)

steel/white steel

commercial grade molded polystyrol electronic control

external digital visual and audible

yes LED light

wire shelves plastic-coated 60 (132.3)

650/530 (25 5/8"/20 7/8") casters with brake at the front, fixed

casters at the rear ergonomic slimline handle


recessed yes

right, reversible on site

50 (2")

113/103 (249.1/227.1)

7113485



 $<sup>^{\</sup>star\star}$  Noise level – sound power re 1pW – according to IEC 60704-2-14.