


Fits in. Stands out.

Take advantage of these special offers!

Packages designed around your applications.
Order your centrifuge, rotor and accessories
with 1 part number.

- Cell Culture
- Microplate and Molecular Biology
- Multi-Application
- Clinical and Healthcare



Thermo Scientific™ Heraeus™ benchtop centrifuges

Fits in.

- To your space, with compact designs
- To your processes and protocols, with complete application-focused packages
- To your future, with a 5-year centrifuge warranty¹ on select packages



Stands out.

Featuring innovative Thermo Scientific™ technologies:

- Secure, push-button rotor exchange in as little as 3 seconds with Auto-Lock™ rotor exchange
- Outstanding ergonomics and performance with Fiberlite™ carbon fiber rotors
- One-handed, certified² sample protection by ClickSeal™ biocontainment lids
- Plus, powerful performance with swing-out capacities up to 4 L and fixed angle rotor speeds over 30,000 x g

¹ Restrictions apply. Please see thermoscientific.com/spinitwarranty for complete warranty offer and coverage details.

² Biocontainment certification by Public Health England, Porton Down, UK.

Cell Culture Packages

Swing-out packages developed for 50 and 15 mL conical tubes. Excellent for cell culture and general life science applications.

Choose for:

- High performance – spin tubes up to 4,816 x g for fast separations
- High capacities from 4 to 40 x 50 mL tubes and 20 to 96 x 15 mL tubes for high throughput
- ClickSeal biocontainment lids designed for user and sample contamination protection
- Optimized deceleration rates for delicate eukaryotic cell pellets, maintaining high yields



Centrifuge	Cat. No.	Package Includes
4 L High Performance and High Throughput		
40 x 50 mL tubes up to 4,200 rpm/4,122 x g (refrigerated) or up to 3,800 rpm/3,374 x g (ventilated)		
Heraeus Multifuge X3FR floor-standing refrigerated centrifuge	75377406	 • TX-1000 rotor with Slide Coat anti-friction finish • TX-1000 aerodynamic buckets (set of 4) • TX-1000 ClickSeal biocontainment lids (set of 4) • Adapters 40 x 50 mL conical • Adapters 96 x 15 mL conical 5-YEAR CENTRIFUGE WARRANTY
Heraeus Multifuge X3F floor-standing ventilated centrifuge	75317406	
Heraeus Multifuge X3R refrigerated centrifuge	75167406	
Heraeus Multifuge X3 ventilated centrifuge	75047406	
3 L High Performance		
28 x 50 mL tubes up to 4,700 rpm/4,816 x g		
Heraeus Multifuge X3R refrigerated centrifuge	75163839	 • TX-750 rotor with Slide Coat anti-friction finish • TX-750 buckets (set of 4) • TX-750 ClickSeal biocontainment lids (set of 4) • Adapters 28 x 50 mL conical • Adapters 56 x 15 mL conical 5-YEAR CENTRIFUGE WARRANTY
Heraeus Multifuge X3 ventilated centrifuge	75013839	
1.6 L High Performance		
16 x 50 mL tubes up to 5,000 rpm/4,696 x g		
Heraeus Multifuge X1R refrigerated centrifuge	75518382	 • TX-400 rotor with Slide Coat anti-friction finish • TX-400 buckets (set of 4) • TX-400 ClickSeal biocontainment lids (set of 4) • Adapters 16 x 50 mL conical • Adapters 36 x 15 mL conical 5-YEAR CENTRIFUGE WARRANTY
Heraeus Multifuge X1 ventilated centrifuge	75118382	
580 mL High Performance		
8 x 50 mL tubes up to 4,500 rpm/3,260 x g		
Heraeus Megafuge 8R refrigerated centrifuge	75130395	 • TX-150 rotor • TX-150 oval buckets (set of 4) • Adapters 8 x 15 mL conical 2-YEAR CENTRIFUGE WARRANTY
Heraeus Megafuge 8 ventilated centrifuge	75100395	

5-YEAR
CENTRIFUGE
WARRANTY¹

5-YEAR
CENTRIFUGE
WARRANTY¹

5-YEAR
CENTRIFUGE
WARRANTY¹

2-YEAR
CENTRIFUGE
WARRANTY²

¹ Restrictions apply. Please see [thermoscientific.com/heraeusspinwarranty](https://www.thermoscientific.com/heraeusspinwarranty) for complete warranty offer and coverage details.

² Subject to Thermo Fisher Scientific's standard limited warranty. See [thermoscientific.com](https://www.thermoscientific.com) or your sales representative for details.

Multi-Application Packages

Multi-rotor packages developed for multi-application environments, meeting your lab's evolving needs with outstanding flexibility and performance.

Choose for:

- Auto-Lock rotor exchange allows you to safely change rotors and applications in as little as 3 seconds with just a push of a button
- High-speed performance for quick spins and fast results
- Large selection of optional rotors and accessories to accommodate a wide variety of tubes
- Fixed angle flexibility with a choice of standard aluminum or high-performance, high-capacity Fiberlite carbon fiber rotors



Centrifuge

Cat. No.

Package Includes

4 L High Throughput Swing-Out and High-Capacity Fixed Angle

40 x 50 mL tubes up to 4,200 rpm/4,122 x g (refrigerated) or 3,800 rpm/3,374 x g (ventilated), and 6 x 250 mL tubes up to 10,000 rpm/15,317 x g

Heraeus Multifuge X3FR floor-standing refrigerated centrifuge 75371762

Heraeus Multifuge X3F floor-standing ventilated centrifuge 75311762

Heraeus Multifuge X3R refrigerated centrifuge 75161762

Heraeus Multifuge X3 ventilated centrifuge 75011762



- TX-1000 rotor with Slide Coat anti-friction finish
- TX-1000 aerodynamic buckets (set of 4)
- TX-1000 ClickSeal biocontainment lids (set of 4)
- Adapters 40 x 50 mL conical
- Adapters 96 x 15 mL conical
- Fiberlite F14-6x250 LE rotor

5-YEAR
CENTRIFUGE
WARRANTY¹

3 L High Performance Swing-Out and Mid-Capacity Fixed Angle

28 x 50 mL tubes up to 4,700 rpm/4,816 x g, and 6 x 100 mL tubes up to 14,500 rpm/24,446 x g

Heraeus Multifuge X3R refrigerated centrifuge 75160063



- TX-750 rotor with Slide Coat anti-friction finish
- TX-750 buckets (set of 4)
- TX-750 ClickSeal biocontainment lids (set of 4)
- Adapters 28 x 50 mL conical
- Adapters 56 x 15 mL conical
- Fiberlite F15-8x50cy rotor

5-YEAR
CENTRIFUGE
WARRANTY¹

1.6 L High Performance and Flexibility

16 x 50 mL tubes up to 5,000 rpm/4,696 x g

Heraeus Multifuge X1R refrigerated centrifuge 75512024

Heraeus Multifuge X1 ventilated centrifuge 75112024



- TX-400 rotor with Slide Coat anti-friction finish
- TX-400 buckets (set of 4)
- TX-400 ClickSeal biocontainment lids (set of 4)
- Adapters 16 x 50 mL conical
- Adapters 36 x 15 mL conical
- HIGHConic II 6 x 100 rotor
- Adapters 6 x 50 mL conical
- M-20 6 x Microplate rotor
- MicroClick 30 x 2 rotor

5-YEAR
CENTRIFUGE
WARRANTY¹

580 mL High Performance and Flexibility

4 x 50 mL tubes up to 4,500 rpm/3,260 x g

Heraeus Megafuge 8R refrigerated centrifuge 75130906

Heraeus Megafuge 8 ventilated centrifuge 75100906



- TX-150 rotor
- TX-150 buckets (set of 4)
- TX-150 ClickSeal biocontainment lids (set of 4)
- Adapters 4 x 50 mL conical
- Adapters 8 x 15 mL conical
- HIGHConic III 6 x 50 rotor
- M10 4 x Microplate rotor
- MicroClick 24 x 2 rotor

2-YEAR
CENTRIFUGE
WARRANTY²

¹ Restrictions apply. Please see thermoscientific.com/heraeusspinitwarranty for complete warranty offer and coverage details.


² Subject to Thermo Fisher Scientific's standard limited warranty. See thermoscientific.com or your sales representative for details.

Microplate and Molecular Biology Packages

Supporting a wide variety of micro-volume solutions with small footprints, high capacities and quiet operation.

Choose for:

- Microplate rotors with capacities from 4 to 24 standard microplates, many with optional certified² biocontainment
- ClickSeal biocontainment lids allow access to samples in seconds
- MicroClick rotors feature extended necks, designed specifically for nucleic acid purification kits
- High speed performance up to 30,279 x g for fast spins and outstanding separations



Centrifuge

Cat. No.

Package Includes

High Throughput Microplate

24 standard microplates up to 4,200 rpm/4,122 x g (refrigerated) or 3,800 rpm/3,374 x g (ventilated)

Heraeus Megafuge 40R refrigerated centrifuge

75191703

Heraeus Megafuge 40 ventilated centrifuge

75041703



- TX-1000 rotor with Slide Coat anti-friction finish
- TX-1000 aerodynamic buckets (set of 4)
- TX-1000 ClickSeal biocontainment lids (set of 4)
- Microplate adapters (set of 4)

5-YEAR
CENTRIFUGE
WARRANTY¹

High Performance Microplate

16 standard microplates up to 4,700 rpm/4,816 x g

Heraeus Multifuge X3R refrigerated centrifuge

75160817

Heraeus Multifuge X3 ventilated centrifuge

75010817



- TX-750 rotor with Slide Coat anti-friction finish
- TX-750 Microplate buckets (set of 4)

5-YEAR
CENTRIFUGE
WARRANTY¹

Standard Microplate

4 standard microplates up to 4,400 rpm/2,576 x g

Heraeus Megafuge 8R refrigerated centrifuge

75130623

Heraeus Megafuge 8 ventilated centrifuge

75100623



- M10 rotor

2-YEAR
CENTRIFUGE
WARRANTY³

1.6 L High Throughput Molecular Biology

16 x 50 mL tubes up to 5,000 rpm/4,696 x g and 30 microtubes up to 14,000 rpm/21,693 x g

Heraeus Multifuge X1R refrigerated centrifuge

75518119

Heraeus Multifuge X1 ventilated centrifuge

75118119



- TX-400 rotor with Slide Coat anti-friction finish
- TX-400 buckets (set of 4)
- TX-400 ClickSeal biocontainment lids (set of 4)
- Adapters 16 x 50 mL conical
- Adapters 36 x 15 mL conical
- MicroClick 30 x 2 rotor (prep kit compatible)

5-YEAR
CENTRIFUGE
WARRANTY¹

580 mL High Performance Molecular Biology

4 x 50 mL tubes up to 4,500 rpm/3,260 x g and 24 microtubes up to 17,850 rpm/30,279 x g

Heraeus Megafuge 8R refrigerated centrifuge

75130115

Heraeus Megafuge 8 ventilated centrifuge

75100115



- TX-150 rotor
- TX-150 buckets (set of 4)
- TX-150 ClickSeal biocontainment lids (set of 4)
- Adapters 4 x 50 mL conical
- Adapters 8 x 15 mL conical
- MicroClick 24 x 2 rotor

2-YEAR
CENTRIFUGE
WARRANTY³

¹ Restrictions apply. Please see thermoscientific.com/heraeusspinitwarranty for complete warranty offer and coverage details.

² Biocontainment certification by Public Health England, Porton Down, UK.

³ Subject to Thermo Fisher Scientific's standard limited warranty. See thermoscientific.com or your sales representative for details.

Clinical and Healthcare Packages

Outstanding reliability with new durable components for super high use laboratories. Ergonomic designs and large user interface allow for ease-of-use.

Choose for:

- High capacity buckets optimized for standard blood tube capacities from 24 to 196 x 5/7 mL tubes
- Auto-Lock rotor exchange for quick rotor removal and easy cleaning
- ClickSeal biocontainment lids help ensure operator safety while spinning patient samples
- Fast acceleration and deceleration times as low as 10 seconds from 0 to 3000 rpm for fast processing and diagnosis



Centrifuge	Cat. No.	Package Includes
4 L High Performance and High Throughput		
196 x 5/7 mL tubes up to 4,200 rpm/4,122 x g (refrigerated) or 3,800 rpm/3,374 x g (ventilated)		
Heraeus Multifuge X3FR floor-standing refrigerated centrifuge	75371771	 
Heraeus Multifuge X3F floor-standing ventilated centrifuge	75311771	
Heraeus Megafuge 40R refrigerated centrifuge	75977271	
Heraeus Megafuge 40 ventilated centrifuge	75477271	
3 L High Throughput		
108 x 5/7 mL tubes up to 4,700 rpm/4,816 x g		
Heraeus Megafuge 40R refrigerated centrifuge	75191923	 
Heraeus Megafuge 40 ventilated centrifuge	75041923	
1.6 L Mid Throughput		
76 x 5/7 mL tubes up to 5,000 rpm/4,696 x g		
Heraeus Megafuge 16R refrigerated centrifuge	75718180	 
Heraeus Megafuge 16 ventilated centrifuge	75318180	
580 mL High Performance		
24 x 5/7 mL tubes up to 4,500 rpm/3,260 x g		
Heraeus Megafuge 8R refrigerated centrifuge	75130102	
Heraeus Megafuge 8 ventilated centrifuge	75100102	
• TX-1000 rotor with Slide Coat anti-friction finish • TX-1000 aerodynamic buckets (set of 4) • TX-1000 ClickSeal biocontainment lids (set of 4) • Adapters 196 x 5/7 mL blood collection • Adapters 148 x 10 mL blood collection 		
• TX-750 rotor with Slide Coat anti-friction finish • TX-750 buckets (set of 4) • TX-750 ClickSeal biocontainment lids (set of 4) • Adapters 108 x 5/7 mL blood collection • Adapters 84 x 10 mL blood collection 		
• TX-400 rotor with Slide Coat anti-friction finish • TX-400 buckets (set of 4) • TX-400 ClickSeal biocontainment lids (set of 4) • Adapters 76 x 5/7 mL blood collection • Adapters 56 x 10 mL blood collection 		
• TX-150 rotor • TX-150 buckets (set of 4) • TX-150 ClickSeal biocontainment lids (set of 4) • Adapters 24 x 5/7 mL blood collection • Adapters 16 x 10 mL blood collection 		

¹ Restrictions apply. Please see thermoscientific.com/heraeusspinitwarranty for complete warranty offer and coverage details.

² Subject to Thermo Fisher Scientific's standard limited warranty. See thermoscientific.com or your sales representative for details.

Increase your package flexibility even further with add-ons

After you choose your application-focused package, simply select from complementary Thermo Scientific products to further enhance your lab's processing capabilities.

Add-on Rotors and Adapters

Thermo Scientific Rotor/Adapter	Cat. No.	Corresponding Centrifuge/Maximum RCF				
		Heraeus Multifuge X3	Heraeus Megafuge 40	Heraeus Multifuge X1	Heraeus Megafuge 16	Heraeus Megafuge 8
Fiberlite F14-6x250 LE rotor	75773662	15,317 x g				
Fiberlite F13-14x50cy rotor	75773661	14,636 x g		12,359 x g		
Fiberlite F15-6x100y rotor	75773698	24,652 x g	18,516 x g	24,652 x g	18,516 x g	
Fiberlite F15-8x50cy rotor	75773663	24,446 x g		22,789 x g		
HIGHConic II rotor	75773620	15,090 x g	15,090 x g	15,090 x g	15,090 x g	
Microliter 48 x 2 rotor	75773602	25,314 x g	25,314 x g	25,314 x g	25,314 x g	
Microliter 30 x 2 rotor	75773652	25,830 x g	25,830 x g	25,830 x g	25,830 x g	
MicroClick 30 x 2 rotor (prep kit compatible)	75775719	21,693 x g	21,693 x g	21,693 x g	21,693 x g	21,693 x g
TX-1000 microplate adapters ⁴	75777303	>3,374 x g	>3,374 x g			
M-20 rotor with standard microplate carriers	75773624	2,272 x g	2,272 x g	2,272 x g	2,272 x g	
Sealed microplate carriers for M-20 rotor ⁴	75773625	2,272 x g	2,272 x g	2,272 x g	2,272 x g	
M10 rotor	75775706					2,576 x g
Standard microplate carriers for M10 rotor ⁴	75775723					2,576 x g
MicroClick 18 x 5 rotor	75775765					22,351 x g

⁴ Rotor sold separately.

Add-on Centrifuges

Thermo Scientific Centrifuge	Includes	Cat. No.	Max Speed	Max RCF
Heraeus Fresco 21 refrigerated microcentrifuge	24 x 1.5/2.0 mL rotor with ClickSeal biocontainment lid	75772426	14,800 rpm	21,100 x g
Heraeus Pico 21 ventilated microcentrifuge	24 x 1.5/2.0 mL rotor with ClickSeal biocontainment lid	75772416	14,800 rpm	21,100 x g
Heraeus Fresco 21 refrigerated microcentrifuge	Dual row 18 x 2.0/0.5 mL rotor with screw-on lid	75772556	14,800 rpm	20,800 x g
Heraeus Pico 21 ventilated microcentrifuge	Dual row 18 x 2.0/0.5 mL rotor with screw-on lid	75772554	14,800 rpm	20,800 x g
mySPIN 12 mini centrifuge	12 x 1.5/2.0 mL tube rotor 32 x 0.2 mL tube rotor (singles or 4 x 8 strip) 0.2 mL and 0.5 mL tube adapters (sets of 12)	75004081	12,500 rpm	9,800 x g
mySPIN 6 mini centrifuge	6 x 1.5/2.0 mL tube rotor 16 x 0.2 mL tube rotor (singles or 2 x 8 strip) 0.2 mL and 0.5 mL tube adapters (sets of 6)	75004061	6,000 rpm	2,000 x g

Valid while supplies last. All promotions valid for orders placed between January 1, 2016 and June 30, 2016 by residents of the United States and Canada. Packages and discounted add-on products must appear on same order to receive discount. "Savings" or "Discount" indicated in any promotion herein represent the maximum percentage difference between the aggregate price for individual purchase of each item contained in the promotional package at the current standard MSRP and the MSRP for the promotional price for each promotional package. Valid only in connection with purchase of promotional packages identified herein using the part numbers listed. All pricing information and values based on current list prices. All offers subject to Thermo Fisher standard terms and conditions of sale, and Thermo Fisher makes no additional warranties, express or implied, in connection with promotional products. Promotional products will be shipped with purchased products. Please see www.thermoscientific.com/spinwarranty for complete warranty offer and coverage details. Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. See representative for additional information. Please note that the value of the specific offer(s) the customer may receive from the manufacturer under this program is a "discount or other reduction in price" to the customer under section 11288(b)(3)(A) of the Social Security Act [42 U.S.C. 1320a-7b (b)(3)(a)]. Accordingly, the customer shall disclose this and any other discounts or other reduction in price received under this program under any state or federal program which provides cost or charge-based reimbursements to the customer for the products or services purchased under this program.



Thermo Scientific™ mySPIN™ 12 and 6 mini centrifuges

thermoscientific.com/spinit

© 2016 Thermo Fisher Scientific Inc. All rights reserved. HERAEUS is a registered trademark of Heraeus Holding GmbH licensed to Thermo Fisher Scientific. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Australia +61 39757 4300
Austria +43 1 801 40 0
Belgium +32 53 73 42 41
China +800 810 5118 or
+400 650 5118
France +33 2 2803 2180
Germany national toll free 0800 1 536 376
Germany international +49 6184 90 6000

India toll free 1800 22 8374
India +91 22 6716 2200
Italy +39 02 95059 552
Japan +81 3 5826 1616
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Nordic/Baltic/CIS countries
+358 9 3291 0200

Russia +7 812 703 42 15
Spain/Portugal +34 93 223 09 18
Switzerland +41 44 454 12 22
UK/Ireland +44 870 609 9203
USA/Canada +1 866 984 3766

Other Asian countries +852 2885 4613
Countries not listed +49 6184 90 6000

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand