

# PET SCALES

**Detecto has them all... large or small!**




**Digital Indicator  
"Locks-On" Weight for Easy Handling**

*Detecto's 758C digital indicator (pictured above on the VET400) is battery powered with automatic sleep mode and automatic activation. This 0.7" high LCD weight display also comes with mounting brackets for desk or wall mounting.*



**Scales for Veterinarians,  
Pet Stores, Kennels or Labs**



# DETECTO PET SCALES

## VET 400

Digital  
Scale up  
to 400 lb



The perfect scale for veterinarians, kennels, labs or anyone handling medium to large size animals. The large, low-profile platform makes it easy for animals to step onto the floor scale without any hesitancy, and an optional rubber mat is also available. The removable stainless steel platform makes clean-up fast and easy. The digital indicator retains the animal's weight even after it is removed from the scale helping ensure simple, accurate operation. The indicator comes with a 8' coiled cable and adjustable mounting bracket for either desk or wall mount convenience.

**Capacity:** 400 lb x .2 lb / 180 kg x .1 kg  
**Platform Size:** 22" W x 42" L x 1 3/4" H  
**Display:** .7" high-contrast LCD  
**Power Source:** 6 "C" size Alkaline batteries or optional AC adapter  
**Construction:** Rugged, steel heavy-duty base, and stainless steel platform  
**Shipping Weight:** 76 lb

## VET 50

Digital  
Scale up  
to 50 lb



This portable, compact-sized scale with a remote weight display is the ideal scale for small animals. Its scale base can be positioned up to 8 feet away from the display. The scale and display are battery powered, and have a power-saving automatic sleep mode. Automatic zero tracking maintains the scale at zero while the hold feature retains the displayed weight after the animal is removed from the scale. The scale's medium base size allows it to be used on the floor or table-top and its light weight makes it ultra portable. An optional rubber mat (as pictured above) is also available.

**Capacity:** 50 lb x .02 lb / 20 kg x .01 kg  
**Platform Size:** 12" W x 14" L x 3" H  
**Display:** .7" high-contrast LCD  
**Power Source:** 6 "C" size Alkaline batteries or optional AC adapter  
**Construction:** Rugged, steel heavy-duty base, and stainless steel platform  
**Shipping Weight:** 15 lb

## 243

Mechanical  
Scale up  
to 20 lb



This sturdy, compact scale incorporates springless balance beam construction. The white baked enamel metal tray holds pets comfortably. The beam is easy to read with blue graduations and numerals on white background. A metric model is also available.

**Capacity:** 20 lb x 2 oz (metric model: 10 kg x 50 g)  
**Tray Size:** 11" W x 19" L x 2 1/2" H  
**Dimensions:** 21" W x 10 1/2" L x 19" H  
**Shipping Weight:** 21 lb

## PS5-2KD


Digital  
Scale up  
to 70 oz



Useful for weighing very small animals, the PS5-2KD features a roomy, stainless steel tray which is removable and easy to clean. Lightweight at just five pounds, the durable PS5-2KD is powered by either four AA batteries (not included) or an optional AC adapter. An easy-to-read 1/2" (13mm) high LCD, and capacity and increments of 70 oz x .1 oz / 1,998 g x 2 g make it ideal for weighing kittens, hamsters, rabbits, lizards, and other small animals.

**Capacity:** 70 oz x .1 oz / 1,998 g x 2 g  
**Tray Size:** 9" L x 6" W x 2" H  
**Dimensions:** 7" L x 7 1/2" W x 2 3/4" D  
**Power Source:** 4 "AA" batteries or optional model with AC adapter and same features  
**Shipping Weight:** 5 lb

Cardinal Scale reserves the right to improve, enhance or modify features and specifications without prior notice. All registered trademarks are the property of their respective owners.



A Division of Cardinal Scale Manufacturing Co.

203 E. Daugherty, Webb City, MO 64870 USA  
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001  
www.detectoscale.com

**SOLD BY:**